

FEMORAL HEAD

Halved or longitudinally halved (bisected) femoral heads are used, for example, in acetabular reconstruction alone or in combination with granules.


Halved femoral heads are available in two different diameters (<45 mm and >45 mm), bisected femoral heads in two different lengths. The height is approx. 20 mm in each case.

Femoral head, halved


Origin: Human
 Tissue: Cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Application: Bone void filler
 Rehydration: Min. 10 minutes

Femoral head, bisected

Origin: Human
 Tissue: Cortico-cancellous bone
 Processing: Allotec® purification procedure
 Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
 Sterilisation: Gamma irradiation
 Application: Bone void filler
 Rehydration: Min. 10 minutes



Bisected Femoral Head – long



Bisected Femoral Head – short

ORDER INFORMATION*


Bisected femoral head

DESCRIPTION	ITEM NUMBER	SIZE
Bisected femoral head	ALO446	short
	ALO447	long



Halved femoral head

DESCRIPTION	ITEM NUMBER	D [mm]	H [mm]
Halved femoral head	ALO441	< 45	20
	ALO444	> 45	20



Halved Femoral Head

* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

BT Medicals AB
 nordics & baltics

Scheelevägen 30-32
 223 63 Lund

+46 364404444
 office@btmedical.com

www.btmedical.com